


Enterostatin (rat)

Chemical Properties

CAS No. :	117137-85-6
Formula:	C ₂₅ H ₄₂ N ₈ O ₈
Molecular Weight:	582.65
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Enterostatin (rat), a procolipase activation peptide with oral activity, selectively diminishes fat consumption and lowers serum cholesterol through a CCK1 receptor-dependent mechanism [1] [2].
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7163 mL	8.5815 mL	17.163 mL
5 mM	0.3433 mL	1.7163 mL	3.4326 mL
10 mM	0.1716 mL	0.8581 mL	1.7163 mL
50 mM	0.0343 mL	0.1716 mL	0.3433 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481