


Substance P-Gly-Lys-Arg

Chemical Properties

CAS No. :	123148-51-6
Formula:	C77H124N24O17S
Molecular Weight:	1690.02
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Substance P-Gly-Lys-Arg, also referred to as β -Preprotachykinin (58-71), is an analog of Substance P.
Targets(IC50)	Neurokinin receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5917 mL	2.9585 mL	5.9171 mL
5 mM	0.1183 mL	0.5917 mL	1.1834 mL
10 mM	0.0592 mL	0.2959 mL	0.5917 mL
50 mM	0.0118 mL	0.0592 mL	0.1183 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481