

Met-RANTES (human)

Chemical Properties

CAS No. :	1883816-50-9
Formula:	C355H543N97O101S6
Molecular Weight:	7972.9
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Met-RANTES (human), a partial antagonist of CCR5, diminishes the infiltration of blood monocytes into the liver [1].
Targets(IC50)	CCR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.1254 mL	0.6271 mL	1.2542 mL
5 mM	0.0251 mL	0.1254 mL	0.2508 mL
10 mM	0.0125 mL	0.0627 mL	0.1254 mL
50 mM	0.0025 mL	0.0125 mL	0.0251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481