

## Eptinezumab

## Chemical Properties

CAS No. : 1644539-04-7

Formula:

Molecular Weight: 143.27 kDa

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Eptinezumab is a humanised monoclonal antibody targeting CGRP (calcitonin gene-related peptide), which prevents CGRP binding to its receptor and may be used for migraine prevention.
Targets(IC50)	CGRP Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.007 mL	0.0349 mL	0.0698 mL
5 mM	0.0014 mL	0.007 mL	0.014 mL
10 mM	0.0007 mL	0.0035 mL	0.007 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481