

Girentuximab

Chemical Properties

CAS No. : 916138-87-9

Formula:

Molecular Weight: 145.55 kDa

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Girentuximab

Biological Description

Description	Girentuximab (G250) is a potent anti-carbonic anhydrase IX (CAIX) monoclonal antibody with anticancer activity for the study of uroepithelial carcinoma PET (ZiPUP).
Targets(IC50)	Carbonic Anhydrase
In vitro	Girentuximab is a chimeric version of a mouse monoclonal antibody (mAb) with IgG1 κ light chain and is capable of recognizing the CAIX antigen [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0069 mL	0.0344 mL	0.0687 mL
5 mM	0.0014 mL	0.0069 mL	0.0137 mL
10 mM	0.0007 mL	0.0034 mL	0.0069 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481