

## Mosunetuzumab

## Chemical Properties

CAS No. : 1905409-39-3

Formula:

Molecular Weight: 146.27 kDa

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Mosunetuzumab (BCT-4465A) is a humanized antibody based on immunoglobulin G1 that targets CD20 (B cells) and CD3 (T cells). Mosunetuzumab redirects T cells to engage and eliminate malignant B cells and can be used to study relapsed or refractory (R/R) B-cell non-Hodgkin lymphoma (B-NHL) and refractory follicular lymphoma.
Targets(IC50)	Others,Immunology/Inflammation related,Integrin

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.0342 mL	0.0684 mL
5 mM	0.0014 mL	0.0068 mL	0.0137 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481