

## Odronextamab

### Chemical Properties

CAS No. : 1801338-64-6

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Odronextamab, a fully human, hinge-stabilized IgG4-based bispecific antibody, targets CD3 receptors on T cells and CD20 receptors on B cells [1].
Targets(IC50)	Immunology/Inflammation related,Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481