

## Utomilumab

## Chemical Properties

CAS No. : 1417318-27-4

Formula:

Molecular Weight: 142.77 kDa

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Utomilumab

## Biological Description

Description	Utomilumab (PF 05082566) is a humanized IgG2 antibody that targets the T-cell co-stimulatory receptor 4-1BB/CD137 and acts as a 4-1BB/CD137 agonist. It exhibits antitumor activity and may be used in the study of patients with advanced malignancies.
Targets(IC50)	TNF

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.007 mL	0.035 mL	0.070 mL
5 mM	0.0014 mL	0.007 mL	0.014 mL
10 mM	0.0007 mL	0.0035 mL	0.007 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481