

Sarilumab

Chemical Properties

CAS No. : 1189541-98-7

Formula:

Molecular Weight: 144.12 kDa

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Sarilumab

Biological Description

Description	Sarilumab (SAR-153191) is a potent interleukin-6 IL-6 receptor antagonist and a human immunoglobulin G1 monoclonal antibody. Sarilumab inhibits IL-6-mediated signaling and reduces inflammation. Sarilumab is used for the prevention and treatment of rheumatoid arthritis.
Targets(IC50)	IL Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0069 mL	0.0347 mL	0.0694 mL
5 mM	0.0014 mL	0.0069 mL	0.0139 mL
10 mM	0.0007 mL	0.0035 mL	0.0069 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481