

## Inotuzumab

## Chemical Properties

CAS No. : 1660159-36-3

Formula:

Molecular Weight: 146.48 kDa

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Inotuzumab (Humanized Anti-CD22 Recombinant Antibody) is a humanized IgG4 antibody targeting CD22. This naked antibody can be used to synthesize antibody-active molecule conjugators [Inotuzumab ozogamicin] and to study acute lymphoblastic leukemia and non-Hodgkin's lymphoma.
Targets(IC50)	Apoptosis,ADC Antibody,Immunology/Inflammation related

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.0341 mL	0.0683 mL
5 mM	0.0014 mL	0.0068 mL	0.0137 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481