

Sibrotuzumab

Chemical Properties

CAS No. :	216669-97-5
Formula:	
Molecular Weight:	148.66 kDa
Storage:	Store at low temperature -20°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Sibrotuzumab (BIBH 1) is a humanized IgG1 monoclonal antibody targeting Fibroblast Activating Protein (FAP) and is used in the study of colorectal cancer and non-small cell lung cancer.
Targets(IC50)	Others,Immunology/Inflammation related

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0067 mL	0.0336 mL	0.0673 mL
5 mM	0.0013 mL	0.0067 mL	0.0135 mL
10 mM	0.0007 mL	0.0034 mL	0.0067 mL
50 mM	0.0001 mL	0.0007 mL	0.0013 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481