

## Galiximab

## Chemical Properties

CAS No. :	357613-77-5
Formula:	
Molecular Weight:	145.24 kDa
Storage:	Store at low temperature -20°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Galiximab (IDEC 114) is a primate-type anti-CD80 immunoglobulin G1 (IgG1)-like monoclonal antibody with antitumor activity that induces in vivo growth inhibition and prolongs survival of B-NHL tumor xenografts. Galiximab has been used in the study of hematologic malignancies.
Targets(IC50)	Apoptosis,Immunology/Inflammation related

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0069 mL	0.0344 mL	0.0689 mL
5 mM	0.0014 mL	0.0069 mL	0.0138 mL
10 mM	0.0007 mL	0.0034 mL	0.0069 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481