

Lokivetmab

Chemical Properties

CAS No. : 1533403-95-0

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lokivetmab, an Anti-Canine IL31 Recombinant Antibody, is a monoclonal antibody targeting canine IL-31, utilized in research related to atopic dermatitis (AD) in dogs [1].
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481