

Pagibaximab

Chemical Properties

CAS No. : 595566-61-3

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pagibaximab is a human chimeric monoclonal antibody targeting lipoteichoic acid, a surface fatty acid of Staphylococcus, for preventing Staphylococcus infections.
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481