

## Quilizumab

## Chemical Properties

CAS No. : 1228538-47-3

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Quilizumab (Anti-Human NGcGM3 Recombinant Antibody) is a humanised monoclonal antibody targeting the M1-prime fragment of membrane-bound IgE, depleting memory B cells and reducing serum IgE levels. It can be used for asthma research.
Targets(IC50)	Apoptosis,Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481