

Vepalimomab

Chemical Properties

CAS No. : 195158-85-1
Formula:
Molecular Weight:
Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Vepalimomab, an Anti-Human AOC3 Recombinant Antibody, is a monoclonal antibody that inhibits Vascular Adhesion Protein 1 (VAP-1), exhibiting anti-inflammatory properties [1].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481