

Etesevimab

Chemical Properties

CAS No. : 2423948-94-9

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Etesevimab (JS016) is a recombinant human IgG1 monoclonal antibody that neutralizes SARS-CoV-2.
Targets(IC50)	Angiotensin-converting Enzyme (ACE),SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481