

Ropeginterferon alfa-2b

Chemical Properties

CAS No. : 1335098-50-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ropeginterferon alfa-2b (Ropeginterferon alfa-2b-njft), a monopegylated interferon-alpha (IFN- α), is used in the study of myeloproliferative neoplasms [1].
Targets(IC50)	Apoptosis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481