

Loncastuximab tesirine

Chemical Properties

CAS No. : 1879918-31-6

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Loncastuximab tesirine, a human cluster of differentiation 19 (CD19)-directed antibody-drug conjugate (ADC), binds to CD19 on the cell membrane, facilitating rapid internalization and inducing cell apoptosis. This mechanism makes it a valuable asset in researching diffuse large B-cell lymphoma [1] [2].
Targets(IC50)	Apoptosis, Antibody-Drug Conjugates (ADCs)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481