

Gresonitamab

Chemical Properties

CAS No. : 2413817-97-5

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Gresonitamab (AMG 910), a half-life extended (HLE) bispecific T-cell engager (BiTE) antibody, specifically targets CD3-positive T cells and CLDN18.2-expressing tumor cells, offering potential research applications in adenocarcinoma [1].
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481