

## Vobarilizumab

## Chemical Properties

CAS No. : 1628814-88-9

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Vobarilizumab (ALX-0061) is a humanized bispecific antibody targeting IL-6R and serum albumin (ALB), and can be used for research on autoimmune diseases such as rheumatoid arthritis.
Targets(IC50)	IL Receptor, Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481