

Sonelokimab

Chemical Properties

CAS No. : 1414386-05-2
Formula:
Molecular Weight:
Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sonelokimab (ALX 0761; M 1095) is a trivalent nanobody combining monovalent camelid-derived nanobodies targeting human interleukin (IL)-17A and IL-17F with human serum albumin VHHs, presenting it as a promising candidate for investigating plaque-type psoriasis treatments [1] [2].
Targets(IC50)	MMP,Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481