

## Abicipar pegol

## Chemical Properties

CAS No. : 1327278-94-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Abicipar pegol (AGN-150998, MP0112), an anti-VEGF DARPIn molecule, is a novel class of small ankyrin repeat domain-containing proteins designed for high specificity and affinity targeting. This compound effectively inhibits angiogenesis and vascular permeability, and has been applied in ocular inflammation-related disease studies through intravitreal injections, significantly reducing mean retinal thickness and leakage area [1] [2].
Targets(IC50)	VEGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481