

Acasunlimab

Chemical Properties

CAS No. : 2253937-12-9

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Acasunlimab (GEN1046), a bispecific antibody (bsAb) designed to target PD-L1 and 4-1BB, functions by enhancing both T-cell and NK-cell response via conditional 4-1BB stimulation and continuously inhibiting the PD-1/PD-L1 axis. This compound is applicable in cancer research [1] [2].
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481