

Alvircept sudotox

Chemical Properties

CAS No. : 137487-62-8
Formula:
Molecular Weight:
Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Alvircept sudotox, a recombinant CD4 fused with exotoxin A from <i>Pneumonas aeruginosa</i> , is utilized in HIV infection research [1].
Targets(IC50)	HIV Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481