

Vibecotamab

Chemical Properties

CAS No. : 2138442-13-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Vibecotamab (XmAb14045) is a potent bispecific antibody that targets CD123 and CD3, promoting T cell-mediated elimination of CD123-expressing cells. It exhibits antitumor properties and is used in studies on acute myeloid leukemia [1].
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481