

## Tuvirumab

## Chemical Properties

CAS No. : 138660-97-6  
Formula:  
Molecular Weight:  
Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Tuvirumab (OST 577; SDZ-OST 577), a human IgG1 subclass monoclonal antibody, targets HBV surface antigen (HBsAg) by binding with high specificity and affinity ( $K=3.6$ nM), holding promise for research into chronic hepatitis B [1] [2].
Targets(IC50)	HBV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481