

Tosatoxumab

Chemical Properties

CAS No. : 1374419-41-6

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tosatoxumab (AR-301; KBSA301) is a fully human immunoglobulin-class monoclonal antibody that targets and neutralises Staphylococcus aureus (S. aureus) α -toxin (α -haemolysin; Hla), thereby serving as an antibiotic-adjuvant therapy for S. aureus ventilator-associated pneumonia (VAP).
Targets(IC50)	Antibacterial, Antibiotic

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481