

## Teropavimab

## Chemical Properties

CAS No. : 2417213-72-8

Formula:

Molecular Weight:

Storage: Store at -80°C  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Teropavimab (3BNC117-LS), an antibody, is utilized in the research of HIV infection [1].
Targets(IC50)	HIV Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481