

## Finotonlimab

## Chemical Properties

CAS No. : 2350298-85-8

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Finotonlimab is a highly specialized recombinant humanized monoclonal antibody of the immunoglobulin G4 isotype target the human programmed cell death protein 1 receptor with exceptional binding affinity and selectivity, restoration of T-lymphocyte mediated anti-tumor immunity by disrupting the inhibitory signaling axis established through interactions with PD-L1 and PD-L2 ligands across various preclinical and clinical experimental models of advanced solid malignancies.
Targets(IC50)	PD-1/PD-L1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481