

Frexalimab

Chemical Properties

CAS No. : 2515463-86-0

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Frexalimab (SAR441344) is a second-generation humanized monoclonal antibody targeting CD40L. It blocks the CD40/CD40L signaling pathway to reduce inflammation and tissue damage, making it applicable for autoimmune diseases such as primary Sjögren's syndrome, systemic lupus erythematosus, and multiple sclerosis.
Targets(IC50)	Immunology/Inflammation related
In vitro	Frexalimab is a second-generation Fc-modified humanized monoclonal antibody against CD40L, whose Fc modification eliminates platelet activation triggered by immune complex stimulation of FcγRIIIa [2].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481