

Gefurulimab

Chemical Properties

CAS No. : 2456407-94-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Gefurulimab (ALXN-1720) is a human-derived bispecific antibody that targets complement C5 and albumin, specifically binding to C5 to inhibit its activation [1].
Targets(IC50)	Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481