

Goflikicept

Chemical Properties

CAS No. : 2416984-26-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Goflikicept (RPH 104), a hybrid protein, selectively targets and neutralizes circulating IL-1 β and IL-1 α . This compound shows promise for the investigation of ST-segment elevation myocardial infarction (STEMI) [1].
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481