

Imvotamab

Chemical Properties

CAS No. : 2573121-53-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Imvotamab (IGM-2323), a CD20xCD3 bispecific IGM antibody featuring a dual-action mechanism, promotes physiological T cell activation to avoid both over-stimulation and immune function down-regulation. This compound is under development for addressing B-cell malignant tumors, including multiple myeloma (MM) and non-Hodgkin's lymphoma (NHL) [1] [2] [3] [4].
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481