

Inezetamab

Chemical Properties

CAS No. : 2568264-76-4

Formula:

Molecular Weight: 195.98 kDa

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Inezetamab is a bispecific monoclonal antibody functioning as an MSLN inhibitor and CD40 agonist, suitable for investigating solid tumours.
Targets(IC50)	TNF,Mesothelin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0051 mL	0.0255 mL	0.051 mL
5 mM	0.001 mL	0.0051 mL	0.0102 mL
10 mM	0.0005 mL	0.0026 mL	0.0051 mL
50 mM	0.0001 mL	0.0005 mL	0.001 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481