

Ipromlimab

Chemical Properties

CAS No. : 2417649-33-1

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ipromlimab is a highly specific humanized monoclonal antibody against human programmed death receptor 1(PD-1). It can block the interaction between PD-1 and its ligand PD-L1/PD-L2 by accurately recognizing and binding the extracellular domain of PD-1, so it can be used as a pharmacological tool to relieve the immunosuppression of T cells, induce the release of anti-tumor cytokines and study the blocking mechanism of immune checkpoints in experimental research.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481