

## Istiratumab

## Chemical Properties

CAS No. : 1509928-04-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Istiratumab (M-6495), a bispecific monoclonal antibody, targets IGF-1R and ErbB3 receptors and facilitates their degradation via the proteasome pathway, offering potential for cancer research [1] [2] [3].
Targets(IC50)	EGFR,IGF-1R

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481