

Izuralimab

Chemical Properties

CAS No. : 2329669-80-7

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Izuralimab is a bispecific IgG1 antibody targeting both the inducible T-cell costimulator (ICOS/CD278) and PD-1 [1].
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481