

Keliximab

Chemical Properties

CAS No. : 174722-30-6
Formula:
Molecular Weight:
Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Keliximab (SB-210396), a chimeric human/macaque IgG1 anti-CD4 monoclonal antibody, demonstrates a binding affinity for soluble CD4 with a K_i value of 1.0 nM. It inhibits T cell proliferation and IL-2 production, making it a valuable agent for cancer research [1] [2].
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481