

## Lecanemab

## Chemical Properties

CAS No. : 1260393-98-3

Formula:

Molecular Weight: 147.06 kDa

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Lecanemab (BAN-2401) is an IgG1 monoclonal antibody targeting soluble aggregated $\beta$ -amyloid (Abeta) with activity against oligomers, protofibrils, and insoluble protofibrils, and may be useful in the study of Alzheimer's disease (AD)-induced mild cognitive impairment and mild AD dementia.
Targets(IC50)	Beta Amyloid

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.034 mL	0.068 mL
5 mM	0.0014 mL	0.0068 mL	0.0136 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481