

## Lenercept

## Chemical Properties

CAS No. : 156679-34-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Lenercept (Ro 45-2081) is a recombinant fusion protein that combines the soluble TNF-receptor (p55) with the Fc portion of human IgG1 [1].
Targets(IC50)	TNF

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481