

Lenvervimab

Chemical Properties

CAS No. : 2055006-79-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lenvervimab (GC1102) is a recombinant human IgG1-type hepatitis B immunoglobulin designed for the research of hepatitis B virus infection [1] [2].
Targets(IC50)	HBV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481