

Lerdelimumab

Chemical Properties

CAS No. : 285985-06-0

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lerdelimumab (CAT-152) is a recombinant human monoclonal antibody of the IgG4 subclass targeting TGF- β 2. It holds potential as an anti-scarring agent in glaucoma research [1] [2].
Targets(IC50)	TGF-beta/Smad

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481