

## Libivirumab

## Chemical Properties

CAS No. : 569658-79-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Libivirumab (17.1.41), a humanized monoclonal antibody targeting hepatitis B virus (HBV), demonstrates neutralizing activity with half maximal inhibitory concentrations (IC50) of 35 ng/mL for HBsAg and 130 ng/mL for HBeAg, respectively [1].
Targets(IC50)	HBV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481