

## Lorigerlimab

## Chemical Properties

CAS No. : 2416595-46-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Lorigerlimab (MGD019) is a bispecific IgG4 dual-affinity re-targeting (DART) antibody that inhibits PD-1 and CTLA-4 to enhance T-cell responses, currently being investigated for metastatic castration-resistant prostate cancer (mCRPC) [1] [2] [3].
Targets(IC50)	PD-1/PD-L1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481