

## Lulizumab pegol

## Chemical Properties

CAS No. : 1421830-13-8

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Lulizumab pegol (BMS-931699) is an anti-CD28 antibody antagonist that inhibits T-cell proliferation, offering potential therapeutic applications for kidney transplantation and autoimmune disease research [1].
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481