

## Manfidokimab

## Chemical Properties

CAS No. : 2428381-56-8

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Manfidokimab (AK120) is a human monoclonal antibody of the IgG4 subclass that specifically targets the interleukin-4 receptor alpha subunit (IL-4R $\alpha$ ). by binding to this shared subunit, Manfidokimab simultaneously blocks the downstream signaling pathways of both interleukin-4 (IL-4) and interleukin-13 (IL-13), and Manfidokimab can be used in studies about moderate to severe atopic dermatitis.
Targets(IC50)	IL Receptor, Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481