

Melredableukin alfa

Chemical Properties

CAS No. : 2056881-92-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Melredableukin alfa (RO-7049665) is a fusion protein consisting of a human IgG1-κ and a mutated human interleukin 2 (IL2 mutein), connected by a peptide linker, and is being researched for potential applications in autoimmune hepatitis and ulcerative colitis [1].
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481