

## Mipeginterferon alfa-2b

## Chemical Properties

CAS No. : 1619970-39-6

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | Mipeginterferon alfa-2b, an analog of interferon alpha-2b (IFNA2B), features an average substitution of 5 among its 11 amino groups (one N-terminal and 10 lysine N6), and exhibits a relative molecular mass of 40 kDa [1]. |
| Targets(IC50) | IFNAR  |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481