

## Modakafusp alfa

## Chemical Properties

CAS No. : 2254522-19-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Modakafusp alfa (TAK-573) is a humanized anti-CD38 IgG4 monoclonal antibody fused with 2 attenuated IFN $\alpha$ 2b molecules, being researched for potential use in treating multiple myeloma [1].
Targets(IC50)	CD38

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481